

# ALFA WASSERMANN VETERINARY REAGENTS BILE ACIDS

Alfa Wassermann offers a reliable, convenient, accurate reagent for your bile acid testing!



## WHY ADD BILE ACIDS TO YOUR TEST MENU?

Total bile acids are metabolized in the liver and serve as a marker for normal and abnormal liver function. Serum bile acids are increased in patients with liver disease.

- 🐾 Smaller kit size helps eliminate waste and maximize reagent usage for patient samples
- 🐾 Reference ranges available for cats and dogs
- 🐾 No open channel bottle requirements
- 🐾 Calibrators included
- 🐾 Controls available for purchase from Alfa Wassermann

## WHAT ARE THE FEATURES OF AWDT BILE ACIDS?

- 🐾 Sample Volume: 4  $\mu$ L
- 🐾 Reagent Volume: R1 270 R2 90
- 🐾 On-Board Stability:  
30 day stability on AWDT chemistry analyzer
- 🐾 Reportable Range: linear from 4.6 to 161.6  $\mu$ mol/L
- 🐾 Analytical Sensitivity:
  - LoD ( $\mu$ mol/L) 4.6
  - LoQ ( $\mu$ mol/L) 4.6



## ORDERING INFORMATION

DESCRIPTION	REORDER #	TESTS/KIT	REFRIGERATED SHIPMENT & STORAGE REQUIRED	KIT SIZE
Bile Acids Reagent (includes Bile Acids Calibrator)	VET101	45	❄️	1 x 5mL, 1 x 13mL

## RELATED ALFA WASSERMANN TESTS

- |           |       |             |       |
|-----------|-------|-------------|-------|
| 🐾 Albumin | 🐾 ALT | 🐾 Bilirubin | 🐾 GGT |
| 🐾 ALP     | 🐾 AST | 🐾 BUN       |       |

For Use on Alfa Wassermann Veterinary Chemistry Systems