

November 2012

## Important Product Information

- **New Lots – Level 1 & Level 2 Chemistry Controls**
  - **Package Insert and Parameter Stickers**

**TO: All Alfa Wassermann Diagnostic Technologies Customers**

Dear Valued Customer:

### Level 1 & Level 2 Chemistry Controls

Please be aware that Alfa Wassermann Diagnostic Technologies will be releasing a new lot of Level 1 and Level 2 chemistry controls shortly. The information for the new lots is as follows:

Product Name	Reorder Number	Lot Number	Expiration Date	Note
Level 1 Chemistry Control	C1-1	99301	2015-05	Please note that ALT (SA1046) is stable for five days.
Level 2 Chemistry Control	C1-2	99302	2015-05	-

Alfa Wassermann recommends performing a parallel study in adherence with good laboratory practices. If you plan to perform this parallel study, we ask that you sequester a portion of your current QC inventory to perform this study.

### Reagent Box Stickers

Additionally, we ask that you note the bright yellow stickers on reagent boxes. These stickers are designed to point out that there have been changes made to package inserts and parameters. Changes on these documents include:

- *Package inserts*: referencing the intended use for the ACE Axcels insert
- *Parameters*: implementation of a low linearity check that only applies ACE Axcels instrument

Should you have any questions, please feel free to contact the **Alfa Wassermann Diagnostic Technologies Customer Solutions Center Toll Free: 866-419-2532**. Thank you for your continued support of AWDT products.

Kind Regards,



Lauren DiPrima  
Marketing Manager

**ALFA WASSERMANN DIAGNOSTIC TECHNOLOGIES, LLC**  
4 HENDERSON DRIVE  
WEST CALDWELL, NJ 07006  
(973) 882-8630  
Fax: (973) 276-0383  
www.AlfaWassermannUS.com  
ACE® and ACE Alera® are registered trademarks of Alfa Wassermann, Inc

